

Textile EEG Ikon

Wearable EEG based on smart textiles for real-world neuroscience research





Ikon

 Self-Managed Technology

Effortlessly conduct EEG research anywhere with self-administered technology that requires no expert supervision.

Comfortable dry textile sensors

The first device on the market with textile sensors—no need for electrolytic substances like gels, saline, or water

• Tested in clinical populations Demonstrated a 90% success rate in clinical

success rate in clinical studies involving over 250 patients.



Wearable textile-EEG device with frontal sensors, that allows brain monitoring for research applications anytime, anywhere by anyone.

> • Medical Standard Designed in accordance with ISO 13485, with medical CE and FDA clearance anticipated for 2025.

- Connectivity and storage Bluetooth Low Energy (BLE) real-time EEG streaming compatible with LSL and other softwares.
- Battery 9+ hours in BLE streaming with local SD.
- Reliable Active shielding with optimized DRL to improve SNR and

Some applications



Ideal for **large-scale data collection** in everyday environments, paving the way for breakthrough insights and faster EEG biomarker discovery.



A solution for **advancing research in precision medicine** and creating customized treatments, enabling tailored interventions based on individual needs.



Conduct research in **real-world scenarios** such as education, business or automotive technology.



Combined with other biosensors, it enables **seamless multimodal research** across various settings, from e-sports and psychology to pharmacological studies.



Technical overview



Layout designed with sensors over prefrontal areas.



Wearable and ultralight (~100g) EEG headset. Quick and easy set up for real-life scenarios.



Reliable textile-EEG monitoring with 24 bits at 256 Hz for 9+ hours. Bluetooth streaming and/or on-board SD storage.



Washable technology that is easy to clean, transport, store, and reuse.

Hardware specifications

Sensors and headset		
EEG channels	Fp1, Fp2, Af7, Af8, A1, A2. Ground in Fpz. Amplifier REF is in FP1; A1 and A2 used for linked ear re-referencing.	
Type of sensors/ electronics	EEG dry sensors with active shielding and optimized DRL	
Head perimeter	52-72 cm	
Wireless Amplifier		
Sampling rate/resolution	256 SPS at 24 bits	
Bandwidth	DC - 40Hz	
Online/real-time impedance check	Yes	
Integrated sensors	Integrated IMU (9 axis): accelerometer, gyroscope and magnetometer	
Input range and noise	±100 mV, < 1 μVRMS (0.5 - 30 Hz) @ 256Hz	
CMRR / Input impedance	> 100 dB @ 50Hx, >50 GΩ	
Data streaming and store		
Data transmission and range	Low energy Bluetooth 4.2 + EDR with 10 meters in direct sight	
Data backup/files	Direct transfer to tablet/phone + data backup on internal SD card.	
Power		
Battery	Rechargable lipo battery. Charging time <3h	
Autonomy	> 9 h	
General		
Weight	~100g	
Warranty	2 years for amplifier	
Certifications	CE (Directive 2014/53) and FCC (part 15)	

Software specifications

Bitbrain software kit (included with equipment)		
Bitbrain real-time SDK	In C/C++ for Windows and Linux.	
Bitbrain data acquisition and visualization suite	Live visualization, streaming or SD recording, data export in CSV, EDF, and raw data visualization.	
Third parties and real-time I/O	LabStreamingLayer LSL compatibility (Matlab, Python, BCl2000, OpenVibe, NeuroPype, etc).	
Third parties data processing	Matlab (EEGLAB, FieldTrip, BCILAB,etc), Python (MNE, etc) and more.	
Bitbrain software platforms (c	untional)	
	puonar)	
SennsLite	Software for data visualization and recording, with large compatibility with real-time I/O and data processing third parties.	
SennsLite Bitbrain Software Development Kit	Software for data visualization and recording, with large compatibility with real-time I/O and data processing third parties. Software kit consisting of different scripts that allow communication and control of the hardware used. It is a starting point for the development of brain-computer interface applications.	

Bundle includes

- EEG headset and amplifier
- Power supply
- Storage base
- Instructions
- Bitbrain Software KitSkin prep wet wipes
- samples

Support

Basic support is included during the lifetime of the product.

- Email support
- Technical assistance
- Knowledge base

Real-world research and applications





Europe

Zaragoza, Spain

Calle. Sta. Teresa de Jesús, 32, 50006 Zaragoza +34 931 444 823

America

New York, United States 228 E 45th Street. Suite 9E

New York, NY 10017



Email info@bitbrain.com

Website www.bitbrain.com